

ENLIGHTENED



Cerumen & Foreign Body Removal from *BIONIX® Medical Technologies*

CPT Codes

69210 – For potential reimbursement of impacted cerumen.

69200 – For potential reimbursement of foreign body removal from the ear.

30300 – For potential reimbursement of foreign body removal from the nose.



The Lighted Articulating Ear Curette™

- Manipulation – Allows the physician to manipulate the curette to work behind cerumen plug or foreign body eliminating the need for the top down sweeping motion that can abrade the ear canal.
- Visualization – A brilliant white light is projected to the tip of the articulating curette illuminating the ear canal. The magnification lens enhances the view of the ear canal to visualize cerumen.
- Safety – Illumination and magnification improves procedural accuracy, avoiding the risk of injury from blind curettage with stainless steel curettes.




The Lighted Ear Curette™ with Magnification

- Saves Time - Visualization of the ear allows you to quickly and efficiently remove cerumen from the ear canal.
- Convenience - Eliminates the need to use a separate light source or magnification device in an attempt to enhance visualization.
- Single Patient Use Tips - Reduces the risk of any cross contamination as well as eliminating the time, labor and costs associated with re-sterilizing steel curettes.



Lighted Forceps for Foreign Body Removal

- Efficacy – Forceps open to accommodate object as large as 5 mm with gripper teeth for a secure hold.
- Simplicity – Simple spring grip design for easy, single handed operation.
- Variety – Ideal for the nose or ear and a variety of foreign bodies.



Product Offerings

The Lighted Ear Curette™ with Magnification


(50 curettes + 1 Light Source + 1 Magnification Lens)

- 2210** - Lighted AngleLoop®
- 2220** - Lighted MicroLoop®
- 2230** - Lighted WaveCurette™
- 2240** - Lighted FlexLoop®
- 2250** - Lighted CeraSpoon®
- 2260** - Lighted InfantScoop®
- 2270** - Lighted VersaLoop®
- 2280** - Lighted CeraPik™
- 2245** - Variety Pack (Includes all tip styles excluding CeraPik)



Lighted Articulating Ear Curette™

- 2511** - Lighted Articulating VersaScoop™
3mm tip (50 Articulating Ear Curette Tips,
1 Light Source, 1 Magnification Lens)



Lighted Forceps for Foreign Body Removal

- 2750** - Box of 10 Forceps, 1 Light Source, 1 Magnification Lens
- 2775** - 10 Replacement Forceps without the Light Source

2299

Lanyard for Lighted Devices




2200

Extra Light Source



2265

Pack of 5 Magnification Lenses



Cerumen Removal Tips & Techniques

View our Cerumen Removal Tips & Techniques Video online at www.BionixMed.com > Videos & Media > Web Videos



1-800-551-6810 | www.BionixMed.com



©2011 Bionix Medical Technologies. Rev. 4/20/11