


Foreign Body Removal

Foreign objects in the nose must be removed because of the very real possibility of aspirating these objects into the lungs. The design of the Single-use Speculum allows for extraction of foreign objects easily from the nasal cavity. A careful examination will reveal the presence of the foreign body, which is then removed by the doctor under direct visualization, spreading the nares with the disposable speculum and using forceps to grasp and remove the object.

Features and Benefits

The Bionix Disposable Speculum provides an easy to use, single-use alternative to conventional metal nasal specula. The Disposable Speculum:

- *Allows for superior view of nasal cavity*
- *One-time usage saves money on sterilization*
- *Smooth, warm tips are more comfortable and less alarming to patient*
- *Comfortable squeeze-grip allows for easy manipulation*
- *Rounded tips protect against irritation*
- *Convenient - always on hand*
- *Facilitates packing or cauterizing of bleeding noses*



Product Numbers

#9878 - Box of 20

#9877 - Box of 48



©2010 Bionix Medical Technologies

5154 Enterprise Blvd. • Toledo, Ohio 43612
Phone: 800.551.6810 • Fax: 800.455.5678

www.BionixMed.com

